

Obtained from LabCorp:

### Sexually Transmitted Diseases (STDs) & NuSwab™ Tests by Nucleic Acid Amplification

For Detection of *C trachomatis*, *N gonorrhoeae*, *T vaginalis*, Herpes Simplex Virus (HSV) 1&2, *Candida* species, and/or Bacterial Vaginosis (BV)

Submit the Gen-Probe APTIMA® vaginal swab transport with orange label (preferred transport) for amplified detection of *C trachomatis*, *N gonorrhoeae*, *T vaginalis*, herpes simplex virus (HSV) 1&2, *Candida* species, and/or bacterial vaginosis (BV) from vaginal fluid specimens. Collect vaginal sample by contacting the swab to the lower third of the vaginal wall and rotating the swab for 10 to 30 seconds to absorb fluid. Immediately place the swab into transport tube and carefully break the swab shaft against the side of the tube. Alternatively, submit the Gen-Probe APTIMA® unisex swab transport with white label for amplified detection of *C trachomatis*, *N gonorrhoeae*, *T vaginalis*, or HSV 1&2 from endocervical or male urethral specimens. A vaginal specimen is required if the unisex swab is used for amplified detection of bacterial vaginosis and/or *Candida* species. All tests can be performed from a single transport device. Do **not** place more than 1 swab in a single transport device. Do **not** use the white-shafted cleaning swab that comes with the unisex swab to collect specimen. Transport at room temperature or refrigerate after collection.




Gen-Probe APTIMA® Vaginal Swab Transport (preferred for NuSwab™ testing)




Gen-Probe APTIMA® Unisex Swab Transport

white swab = cleaning then discard



blue swab = specimen collection

Submit the Gen-Probe APTIMA® urine transport with yellow label for amplified detection of *C trachomatis*, *N gonorrhoeae*, and/or *T vaginalis* from male or female first-void urine specimens. Total volume of urine submitted must fall between the two black lines indicated on the transport tube label. Transport at room temperature or refrigerate after collection.



Gen-Probe APTIMA® Urine Transport (for Ct, Ng, and Tv testing only)

