


Hunterdon Medical Center Laboratory

Order: 100997

Respiratory 4 Viral Panel by PCR: Collection Kit

	Nasopharyngeal Swab and Transport Medium
Catalog #	SWAB/B-100 (Copan UTM® P/N 330C and Copan Nylon Flocked Swab® P/N 503CS01)
Intended Use	Designed to collect, preserve and transport respiratory virus specimen. (Xpert Viral Transport Medium® or UTM-RT®)
Kit Contents	1 flocked NP collection swab 1 tube swab transport medium (3mL, red cap)

<p>Specimens should be transported at 2-8°C</p> <p>Stability: 48 hours Ambient 15-30°C 7 Days Ref 2-8°C</p>	
---	--